

Product Datasheet

Human CD25 ELISA Kit 950.500

Product Name	Human CD25 ELISA Kit
Catalog Number	950.500
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	68.75-2200 pg/ml
Sensitivity	8 pg/ml
Tested Applications	ELISA
Alternative Names	IL-2R
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble receptors
Specificity	Recognizes both natural and recombinant human soluble IL-2R
Sample Types	Serum, Plasma, Cell culture supernatant