Medix Biochemica

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 Fl-02180 Espoo

Finland

Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

Product Datasheet

Human CD25 ELISA Kit 950.500

Product Name Human CD25 ELISA Kit

Catalog Number 950.500

Description

Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed

procedure.

Range 68.75-2200 pg/ml

Sensitivity 8 pg/ml

Alternative Names IL-2R

Brand Diaclone

Storage Store at +4°C

Cross Reactivity No cross reactivity with other human soluble receptors

Specificity Recognizes both natural and recombinant human soluble IL-2R

Sample Types Serum, Plasma, Cell culture supernatant