

Product Datasheet

Human CD138 ELISA Kit 950.640

Product Name	Human CD138 ELISA Kit
Catalog Number	950.640
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	8-256 ng/ml
Sensitivity	4.94 ng/ml
Tested Applications	ELISA
Alternative Names	Syndecan-1
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes natural human soluble syndecan-1
Sample Types	Serum, Plasma, Cell culture supernatant