

Product Datasheet

Anti-h LH 5301 SP-5 100016

Product Name	Anti-h LH 5301 SP-5
Catalog Number	100016
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components
Species/Host	Mouse
Tested Applications	CLIA, LF, FIA
Alternative Names	LH, Lutropin, Interstitial Cell-Stimulating Hormone
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8 °C
Note	Nilsson et al. (2001) analyzed epitopes of 30 different LH mAbs. Antibody 5301 was classified as belonging to epitope group beta 2, recognizing intact LH, its beta subunit as well as a common variant of LH. Antibody 5301 did not cross react with TSH or FSH but did so with hCG.
Isotype	IgG1
Clonality	Monoclonal
Epitope	Beta 2 as described in Nilsson et al. (2001) Two antibodies binding to the same, or closely located epitopes, belong to the same group and hence cannot be used as a pair in a sandwich assay. Epitope group numbering does not give any detailed information where the epitope is located.
Purity	≥ 95 %
Affinity constant	$K_A = 4.4 \times 10^{11} \text{ l/M}$; $K_D = 2.3 \times 10^{-12} \text{ M}$ (= 2.3 pM)
Buffer	0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.3-7.2
Cross Reactivity	LH α 13 %, LH β 170 %, FSH 5 %, hCG 138 %, TSH 0.03 %
Specificity	Antibody recognizes human luteinizing hormone (lutropin), and its beta-subunit
Shelf Life	36 months

- References**
- Federici, M.M., Fraser, R., Lundqvist, C., and Lankford, J.C., (1982) Production and characterization of monoclonal antibodies human lutenizing hormones. Fed. Proc., 41
- Nilsson, C., Seppälä, M., and Pettersson, K., (21) Immunological characterization of human luteinizing hormone with special regard to a common genetic variant. J.Endocrinol. 168:1-116
- Pettersson, K.S.I., and Söderholm J.R-M., (199) Ultrasensitive two-site immunometric assay of human lutropin by time-resolved fluorometry. Clin. Chem. 36(11):1928-1933
- Pettersson, K.S.I., and Söderholm J.R-M., (1991) Individual differences in lutropin immunoreactivity revealed by monoclonal antibodies. Clin. Chem. 37(3):333-34
- Vilja, P., Wichmann, L., Isola, J., and Tuohimaa, P., (1988) Monoclonal-antibody based noncompetitive avidin-biotin assay for lutropin in urine.Clin. Chem. 34(8):1585-159