

Product Datasheet

Anti-h Inhibin alpha 9401 SPTN-5 100711

Product Name	Anti-h Inhibin alpha 9401 SPTN-5
Catalog Number	100711
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Species/Host	Mouse
Tested Applications	FIA
Alternative Names	Inhibin-A
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8 °C
Isotype	Mouse IgG1
Clonality	Monoclonal
Epitope	Amino acids 233-264 (STPLMSWPWSPSALRLLQRPPEEPAAHANCHR)
Purity	≥ 95 %
Affinity constant	N/D
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.3-7.2
Cross Reactivity	Not determined
Specificity	Antibody recognizes human inhibin alpha subunit
Shelf Life	36 months