

Product Datasheet

Anti-C. difficile GDH 11001 SPTN-5 100714

Product Name	Anti-C. difficile GDH 11001 SPTN-5
Catalog Number	100714
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	N/D
Associated Products	Recombinant C. difficile GDH antigen, Medix Biochemica 610121
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	7.0-7.8
Cross Reactivity	N/D
Specificity	Antibody recognizes Clostridium difficile glutamate dehydrogenase
Matched Pairs	Forms pair with itself
Shelf Life	18 months