

Product Datasheet

Aprotinin, Bovine Lung 125-10-1,000

Product Name	Aprotinin, Bovine Lung
Catalog Number	125-10-1,000
Alternative Names	Pancreatic Trypsin Inhibitor, Basic Protease Inhibitor, BPI, BPTI
Brand	Lee Biosolutions
Form/Appearance	Lyophilized White to off-white hygroscopic powder
Concentration	≥ 0.9 mg protein/mg solid
Storage	2-8° C
Purity	≥ 95% (HPLC)
CAS Number	9087-70-1
Source	Bovine Lung
Molecular Weight	6511.44
Shelf Life	2 years
Biological Activity	3.0 - 8.0 TIU/mg solid