

Product Datasheet

Anti-h NT-proBNP 1312 SPTN-5 100717

Product Name	Anti-h NT-proBNP 1312 SPTN-5
Catalog Number	100717
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Species/Host	Mouse
Tested Applications	LF, CLIA, FIA
Alternative Names	N-Terminal Prohormone of Brain Natriuretic Peptide, BNPT
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Amino acids 59-70 (EGIRGHRKMVLY)
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.6-7.3
Cross Reactivity	Does not recognize canine NT-proBNP
Specificity	Antibody recognizes human N-terminal pro-B-type natriuretic peptide
Matched Pairs	Capture - detection CLIA: 1306 - 1309, 1306 - 1311, 1306 - 1312, 1308 - 1309, 1308 - 1310, 1308 - 1311 and 1308 - 1312 LF: 1306 - 1309, 1306 - 1310 FIA: 1306 - 1310, 1306 - 1312, 1309 - 1307, 1310 - HM210, HM071 - 1309, 1310 - HM071, 1307 - 1311
Shelf Life	Unspecified