

Product Datasheet

Anti-h Galectin-3 10302 SPTN-5 100731

Product Name	Anti-h Galectin-3 10302 SPTN-5
Catalog Number	100731
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	CLIA, FIA
Alternative Names	Mac-2 antigen, IgE-binding protein, L-29, CBP30
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Located in the amino acid region A66 - T104
Purity	≥ 95 %
Associated Products	Recombinant Galectin-3 antigen, Medix Biochemica 610144
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.1-6.6
Cross Reactivity	Does not recognize Galectin-1, Galectin-9, or Galectin-3-binding protein
Specificity	Antibody recognizes human galectin-3
Shelf Life	Unspecified