

Product Datasheet

Anti-h Galectin-3 10303 SPTN-5
100732

| | |
|---------------------|--|
| Product Name | Anti-h Galectin-3 10303 SPTN-5 |
| Catalog Number | 100732 |
| Description | Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components. |
| Tested Applications | FIA, CLIA |
| Alternative Names | Mac-2 antigen, IgE-binding protein, L-29, CBP30 |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid, may turn slightly opaque during storage |
| Concentration | 5.0 mg/ml (+/- 10 %) |
| Storage | +2-8°C |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Epitope | Located in the amino acid region A66 - T104 |
| Purity | ≥ 95 % |
| Associated Products | Recombinant Galectin-3 antigen, Medix Biochemica 610144 |
| Buffer | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative |
| IEF Profile | 6.5-7.1 |
| Cross Reactivity | Does not recognize Galectin-1, Galectin-9, or Galectin-3-binding protein |
| Specificity | Antibody recognizes human galectin-3 |
| Shelf Life | Unspecified |