

Product Datasheet

Anti-h Lp-PLA2 11207 SPTN-5 100753

Product Name	Anti-h Lp-PLA2 11207 SPTN-5
Catalog Number	100753
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Species/Host	Mouse
Tested Applications	FIA
Alternative Names	Platelet-activating factor acetylhydrolase (PAF-AH)
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8 °C
Note	Casein is recommended to be used in blocking buffer instead of BSA.
Isotype	IgG1
Clonality	Monoclonal
Epitope	Not determined
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.0-6.8
Cross Reactivity	Not determined
Specificity	Antibody recognizes human lipoprotein-associated phospholipase A ₂
Shelf Life	18 months