

Product Datasheet

Anti-h AMH 11303 SPTN-5 100758

Product Name	Anti-h AMH 11303 SPTN-5
Catalog Number	100758
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Species/Host	Mouse
Tested Applications	CLIA, LF, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N-terminal region of AMH within amino acids Arg26-Arg451
Purity	≥ 95 %
Affinity constant	N/D
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.4-7.0
Cross Reactivity	Not determined
Specificity	Antibody recognizes human anti-Müllerian hormone (AMH)
Matched Pairs	FIA (capture – detection): 11301 – 11303, 11303 – 11301 and 11309 – 11302 LF (membrane – particles): 11304 – 11303 CLIA (capture – detection): 11303 – 11309, 11304 – 11309 and 11303 – 11302
Shelf Life	36 months