

Product Datasheet

Anti-h AMH 11305 SPTN-5
100760

Product Name	Anti-h AMH 11305 SPTN-5
Catalog Number	100760
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	LF, CLIA, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N-terminal region of AMH within amino acids Arg26-Arg451
Purity	≥ 95 %
Affinity constant	N/D
Associated Products	Recombinant Anti-Mullerian Hormone (AMH) Antigen LA578 Recombinant Anti-Mullerian Hormone (AMH) Antigen LA467
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
IEF Profile	6.6–7.0
Cross Reactivity	N/D
Specificity	Antibody recognizes human anti-Müllerian hormone (AMH)
Matched Pairs	FIA: 11301 (capture) – 11303 (detection), 11303 – 11301 and 11309 – 11302 LF: 11304 (membrane) – 11303 (particles) CLIA: 11303 (capture) – 11309 (detection) and 11304 – 11309
Shelf Life	Unspecified