

Product Datasheet

Anti-h PCT 4008 SPTN-5  
100769

Product Name	Anti-h PCT 4008 SPTN-5
Catalog Number	100769
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	LF, FIA, CLIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG3
Clonality	Monoclonal
Epitope	Anti-h PCT 4008 antibody recognizes the katalcalcin domain of procalcitonin, i.e. the region Asp(96) - Asn(116), as calculated from the peptide sequence of the mature PCT protein, without the signal sequence. The binding epitope is located in the region Asp(96)-Ser(110).
Purity	≥ 95 %
Affinity constant	N/D
Associated Products	Recombinant PCT antigen, Medix Biochemica 610080
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
IEF Profile	8.2-9.1
Cross Reactivity	N/D
Specificity	Antibody recognizes human procalcitonin (PCT)
Shelf Life	Unspecified