

Product Datasheet

Anti-h SAA 2205 SPTN-5 100802

Product Name	Anti-h SAA 2205 SPTN-5
Catalog Number	100802
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	LF, IT, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Recognizes aa 64-71 but also reacts with aa 28-38 indicating that the epitope may be partly conformational.
Purity	≥ 95 %
Associated Products	Native human SAA, EastCoast Bio cat LA449, recombinant human SAA, Medix Biochemica cat 610070 and Lee Biosolutions cat 527-18
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.7 – 7.2
Cross Reactivity	Recognizes human SAA1. Does not recognize human SAA2, human SAA4, feline, canine, equine or bovine SAA.
Specificity	Antibody recognizes human serum amyloid A (SAA)
Shelf Life	Unspecified