

Product Datasheet

Anti-h SAA 2213 SPTN-5 100810

Product Name	Anti-h SAA 2213 SPTN-5
Catalog Number	100810
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Species/Host	Mouse
Tested Applications	LF, FIA, IT
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8 °C
Isotype	IgG2a
Clonality	Monoclonal
Epitope	Conformational
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	7.3 - 8.7
Cross Reactivity	Recognizes human SAA1, human SAA2, feline, canine and equine SAA. Does not recognize human SAA4 or bovine SAA.
Specificity	Antibody recognizes human serum amyloid A (SAA)
Shelf Life	Unspecified