

Product Datasheet

Anti-h Prealbumin 11604 SPTN-5 100831

| | |
|----------------------------|----------------------------------------------------------------------------------------------------|
| Product Name | Anti-h Prealbumin 11604 SPTN-5 |
| Catalog Number | 100831 |
| Description | Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components. |
| Species/Host | Mouse |
| Tested Applications | FIA, IT |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid, may turn slightly opaque during storage |
| Concentration | 5.0 mg/ml (+/- 10 %) |
| Storage | +2-8 °C |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Epitope | Not determined |
| Purity | ≥ 95 % |
| Buffer | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative |
| IEF Profile | 6.7-7.3 |
| Cross Reactivity | Not determined |
| Specificity | Antibody recognizes human prealbumin |
| Shelf Life | Unspecified |