

Product Datasheet

Anti-h GDF-15 4905 SPTN-5
100838

Product Name	Anti-h GDF-15 4905 SPTN-5
Catalog Number	100838
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Species/Host	Mouse
Tested Applications	CLIA, FIA
Alternative Names	Growth/differentiation factor 15, GDF15, Macrophage inhibitory cytokine 1, MIC1, MIC-1
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Not determined
Purity	≥ 95 %
Affinity constant	Not Applicable (N/A)
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
IEF Profile	6.0-7.2
Cross Reactivity	Does not recognize GDF-9, GDF-11, TNF-beta, TGF-beta or CRP
Specificity	Antibody recognizes human growth differentiation factor 15 (GDF-15)
Shelf Life	Unspecified