

Product Datasheet

Propoxyphene (PPX) Antibody 100895

Product Name	Propoxyphene (PPX) Antibody
Catalog Number	100895
Description	Protein A purified monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	≥ 1.5 mg/ml
Storage	+2-8 °C
Note	This product was previously available with Catalog codes P001M and CB302. It has been produced using the same original mouse monoclonal hybridoma cell line.
Isotype	IgG1
Clonality	Monoclonal
Epitope	Not determined
Purity	≥ 95 %
Affinity constant	N/D
Buffer	0.05 M Potassium Phosphate Buffer, pH 7.2 with 0.85 % NaCl, 0.05 % NaN ₃ as a preservative
IEF Profile	6.7 - 7.6
Cross Reactivity	Not determined
Specificity	Antibody recognizes Propoxyphene
Matched Pairs	170032 (BSA conjugate as capture)
Shelf Life	36 months