

Product Datasheet

Propoxyphene (PPX) Antibody 100895

| | |
|----------------------------|--|
| Product Name | Propoxyphene (PPX) Antibody |
| Catalog Number | 100895 |
| Description | Protein A purified monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components. |
| Species/Host | Mouse |
| Tested Applications | ELISA, LF |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid, may turn slightly opaque during storage |
| Concentration | ≥ 1.5 mg/ml |
| Storage | +2-8 °C |
| Note | This product was previously available with Catalog codes P001M and CB302. It has been produced using the same original mouse monoclonal hybridoma cell line. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Epitope | Not determined |
| Purity | ≥ 95 % |
| Affinity constant | N/D |
| Buffer | 0.05 M Potassium Phosphate Buffer, pH 7.2 with 0.85 % NaCl, 0.05 % NaN ₃ as a preservative |
| IEF Profile | 6.7 - 7.6 |
| Cross Reactivity | Not determined |
| Specificity | Antibody recognizes Propoxyphene |
| Matched Pairs | 170032 (BSA conjugate as capture) |
| Shelf Life | 36 months |