

Product Datasheet

Tricyclic Antidepressant (TCA) Antibody 100896

Product Name	Tricyclic Antidepressant (TCA) Antibody
Catalog Number	100896
Description	Protein A purified monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	LF, ELISA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	≥ 1.5 mg/ml
Storage	+2-8°C
Note	This product was previously available with Catalog codes T002M and CB300. It has been produced using the same original mouse monoclonal hybridoma cell line.
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	N/D
Associated Products	T001A or CB500 (BSA conjugate as capture)
Buffer	0.05M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl, 0.05% NaN ₃ as a preservative
IEF Profile	6.1-7.1
Cross Reactivity	N/D
Matched Pairs	T001A or CB500 (BSA conjugate as capture)
Shelf Life	18 months