Medix Biochemica

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 FI-02180 Espoo

Finland

Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

Product Datasheet

Phencyclidine (PCP) Antibody 100900

Product Name

Phencyclidine (PCP) Antibody

Catalog Number 100900

Description

Protein A purified monoclonal mouse antibody, cultured in vitro

under conditions free from animal-derived components.

Species/Host Mouse

Tested **Applications** LF, ELISA

Immunogen

Phencyclidine-BSA

Brand

Medix Biochemica

Form/AppearaLiquid, may turn slightly opaque during storage

Concentratior ≥ 1.5 mg/ml

Storage +2-8 °C

Note

This product was previously available with Catalog code P01-99-64M-P. It has been produced using the same original

mouse monoclonal hybridoma cell line.

lgG1 Isotype

Clonality Monoclonal

Epitope Not determined

≥ 95 % **Purity**

Affinity constant N/D

Buffer 0.05 M Potassium Phosphate Buffer, pH 7.2 with 0.85 % NaCl,

0.05 % NaN3 as a preservative

Cross Reactivity Phenyl cyclohexyl morpholine 4.3%, 4-Hydroxy phencyclidine 1.8%. No reactivity with benzoylecgonine, opiates, cannabinoids,

amphetamines.

Specificity Antibody primarily recognizes phencyclidine

Shelf Life 24 months