

## Product Datasheet

### Tricyclic Antidepressant (TCA)-BSA Antigen 170030

|                        |   |
|------------------------|---|
| <b>Product Name</b>    | Tricyclic Antidepressant (TCA)-BSA Antigen  |
| <b>Catalog Number</b>  | 170030  |
| <b>Description</b>     | Synthetic Tricyclic Antidepressant, BSA conjugated  |
| <b>Brand</b>           | Medix Biochemica  |
| <b>Form/Appearance</b> | Liquid  |
| <b>Concentration</b>   | Lot-specific  |
| <b>Storage</b>         | Short term storage ( $\leq 2$ weeks): 2-8 °C . Long term: -20 °C. Avoid repeated freezing and thawing.                              |
| <b>Note</b>            | This product was previously available with Catalog codes T001A and CB500. It has been produced using the same conjugation protocol. |
| <b>Buffer</b>          | Distilled water   |
| <b>Matched Pairs</b>   | 100896 (detection)  |
| <b>Shelf Life</b>      | 60 months   |