

Product Datasheet

Tricyclic Antidepressant (TCA)-BSA Antigen 170030

Product Name	Tricyclic Antidepressant (TCA)-BSA Antigen
Catalog Number	170030
Description	Synthetic Tricyclic Antidepressant, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	Lot-specific
Storage	Short term storage (≤ 2 weeks): 2–8 °C . Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with Catalog codes T001A and CB500. It has been produced using the same conjugation protocol.
Buffer	Distilled water
Matched Pairs	100896 (detection)
Shelf Life	60 months