

# Product Datasheet

## Propoxyphene (PPX)-BSA Antigen 170032

Product Name	Propoxyphene (PPX)-BSA Antigen
Catalog Number	170032
Description	Synthetic Propoxyphene, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	Lot-specific
Storage	Short term storage ( $\leq$ 2 weeks): 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog codes P001A and CB501. It has been produced using the same conjugation protocol.
Associated Products	100895 (detection)
Buffer	Distilled water
Shelf Life	60 months