

## Product Datasheet

### Morphine-BSA Antigen 170036

<b>Product Name</b>	Morphine-BSA Antigen
<b>Catalog Number</b>	170036
<b>Description</b>	Synthetic Morphine, BSA conjugated
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	Lot-specific
<b>Storage</b>	Short Term ( $\leq 2$ weeks): +2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
<b>Note</b>	This product was previously available with a Catalog code P81-99-01A. It has been produced using the same conjugation protocol.
<b>Associated Products</b>	100897 (detection)
<b>Buffer</b>	Distilled Water
<b>Shelf Life</b>	60 months