

Product Datasheet

Morphine-BSA Antigen 170036

| | |
|---------------------|---|
| Product Name | Morphine-BSA Antigen |
| Catalog Number | 170036 |
| Description | Synthetic Morphine, BSA conjugated |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid |
| Concentration | Lot-specific |
| Storage | Short term storage (≤ 2 weeks): 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing. |
| Note | This product was previously available with a Catalog code P81-99-01A. It has been produced using the same conjugation protocol. |
| Associated Products | 100897 (detection) |
| Buffer | Distilled water |
| Shelf Life | 60 months |