

Product Datasheet

Ketamine-BSA Antigen 170037

Product Name	Ketamine-BSA Antigen
Catalog Number	170037
Description	Synthetic Ketamine, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	Lot-specific
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog code P82-99-31A. It has been produced using the same conjugation protocol.
Buffer	Distilled Water
Matched Pairs	100902 or PP2111F (detection)
Shelf Life	5 years from manufacturing at -20 °C, subject to recertification thereafter