

# Product Datasheet

## Tramadol-BSA Antigen 170038

Product Name	Tramadol-BSA Antigen
Catalog Number	170038
Description	Synthetic Tramadol, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	Lot-specific
Storage	Short term storage ( $\leq 2$ weeks): 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog code P82-99-61A. It has been produced using the same conjugation protocol.
Associated Products	100901 (detection)
Buffer	Distilled water
Shelf Life	60 months