

Product Datasheet

Urobilinogen, Synthetic (Reagent for Urinalysis) 651-10-0.1

Product Name	Urobilinogen, Synthetic (Reagent for Urinalysis)
Catalog Number	651-10-0.1
Alternative Names	Stercobilinogen
Brand	Lee Biosolutions
Form/Appearance	Liquid Amber to orange liquid
Concentration	25 g/L (UBG assay)
Storage	-20 °C
CAS Number	14684-37-8
Source	Synthetic
Molecular Weight	590.72 - 598.78 g/mol
Shelf Life	1 year