

# Product Datasheet

## Methylphenidate-BSA Antigen 170039

|                     |   |
|---------------------|---|
| Product Name        | Methylphenidate-BSA Antigen   |
| Catalog Number      | 170039  |
| Description         | Synthetic Methylphenidate, BSA conjugated   |
| Brand               | Medix Biochemica  |
| Form/Appearance     | Liquid  |
| Concentration       | Lot-specific  |
| Storage             | Short term storage ( $\leq 2$ weeks): 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.                           |
| Note                | This product was previously available with a Catalog code P83-99-01A. It has been produced using the same conjugation protocol. |
| Associated Products | P03-99-01R-IF, P03-99-02R-IF, P03-99-03R-IF (detection)   |
| Buffer              | Distilled water   |
| Shelf Life          | 60 months   |