

Product Datasheet

Methylphenidate-BSA Antigen 170039

Product Name	Methylphenidate-BSA Antigen
Catalog Number	170039
Description	Synthetic Methylphenidate, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	Lot-specific
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog code P83-99-01A. It has been produced using the same conjugation protocol.
Associated Products	P03-99-01R-IF, P03-99-02R-IF, P03-99-03R-IF (detection)
Buffer	Distilled Water
Shelf Life	5 years from manufacturing at -20 °C, subject to recertification thereafter.