

Product Datasheet

Fentanyl-BSA Antigen 170040

Product Name	Fentanyl-BSA Antigen
Catalog Number	170040
Description	Synthetic Fentanyl, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	Lot-specific
Storage	Short term storage (≤ 2 weeks): 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog code P81-99-51A. It has been produced using the same conjugation protocol.
Buffer	Distilled water
Matched Pairs	100998, 100974, P01-99-53R-IF, P01-99-54R-IF, P01-99-55R-IF (detection)
Shelf Life	60 months