

Product Datasheet

Recombinant Human Interleukin-6 (IL-6)-BSA 715-H01-025-BSA

Product Name	Recombinant Human Interleukin-6 (IL-6)-BSA
Catalog Number	715-H01-025-BSA
Description	Recombinant Protein
Species/Host	HEK cells
Alternative Names	NULL
Brand	Diaclone
Form/Appearance	Lyophilized from PBS pH 7.4 (NaCl 155mM, Na ₂ HPO ₄ 8mM & KH ₂ PO ₄ , 1.8mM, Trehalose) + 1% BSA
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles. Stable one week at 2 to 8°C after reconstitution under sterile conditions. Stable 12 months -20°C to -80°C after reconstitution under sterile conditions
Purity	>95% by SDS page
Biological Activity	Measured in a cell based assay on stimulated XG1 cells (EC50= 0.4 ng/ml)
Sequence	Ala28-Met212
Tag	His Tag C-Terminus