

Product Datasheet

Zolpidem-BSA Antigen 170041

Product Name	Zolpidem-BSA Antigen
Catalog Number	170041
Description	Synthetic Zolpidem, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	Lot-specific
Storage	Short term storage (≤ 2 weeks): 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog code P83-99-11A. It has been produced using the same conjugation protocol.
Buffer	Distilled water
Matched Pairs	JP046, JP047, JP048, HM311 or HM312 (detection)
Shelf Life	60 months