

Product Datasheet

Carisoprodol (Soma)-BSA Antigen 170046

Product Name	Carisoprodol (Soma)-BSA Antigen
Catalog Number	170046
Description	Synthetic Carisoprodol, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	Lot-specific
Storage	Short term storage (≤ 2 weeks): 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog code P82-99-71A. It has been produced using the same conjugation protocol.
Buffer	Distilled water
Matched Pairs	P02-99-72R-IF or P02-99-73R-IF (detection)
Shelf Life	60 months