

Product Datasheet

Human IL-2 ELISpot Kit 856.001PC

Product Name	Human IL-2 ELISpot Kit
Catalog Number	856.001PC
Description	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IL-2.
Specificity	Recognizes natural human IL-2
References	Junwei, W. et al., Cancer Cell Int.,2016; 17: 10. - Pubmed link