

Product Datasheet

Human IL-2 ELISpot Set 856.001.001S

Product Name	Human IL-2 ELISpot Set
Catalog Number	856.001.001S
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 1 sterile plate.
Tested Applications	ELISpot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IL-2.
Specificity	Recognizes natural human IL-2