
Product Datasheet

Human IL-4 ELISpot Set 856.011.005S

| | |
|----------------------------|--|
| Product Name | Human IL-4 ELISpot Set |
| Catalog Number | 856.011.005S |
| Description | Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 5 sterile plates. |
| Tested Applications | ELISpot |
| Alternative Names | NULL |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Note | Spots read visually or with reader |
| Cross Reactivity | No cross reactivity with other human cytokines |
| Specificity | Recognizes natural human IL-4 |