

## Product Datasheet

### NGAL Antigen, Recombinant 610012

<b>Product Name</b>	NGAL Antigen, Recombinant
<b>Catalog Number</b>	610012
<b>Description</b>	Recombinant human neutrophil gelatinase-associated lipocalin (NGAL) protein with a C-terminal histidine tag.
<b>Species/Host</b>	E. coli
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	+2-8 °C
<b>Associated Products</b>	Anti-h NGAL 4202, 4203, 4204 and 4205
<b>Buffer</b>	Lyophilized from phosphate buffered saline, pH 7.4; containing 3 % sucrose, 2 % D-mannitol and 0.01 % Tween 20 as stabilizers
<b>Molecular Weight</b>	Predicted molecular weight: 22 kDa
<b>Shelf Life</b>	Recertification 30 months
<b>Expression System</b>	E. coli