
Product Datasheet

Human IL-10 ELISpot Set 856.031

Product Name	Human IL-10 ELISpot Set
Catalog Number	856.031
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 1 non-sterile plate. Spots read visually or with reader.
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10.
Specificity	Recognizes natural human IL-10