

Product Datasheet

Recombinant hCoV OC43 N, 100 µg 610040

Product Name	Recombinant hCoV OC43 N, 100 µg
Catalog Number	610040
Description	Recombinant human coronavirus OC43 nucleoprotein with an N-terminal histidine tag.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Storage	Storage at 2-8 °C. After reconstitution stable for 1 week at 2-8 °C and 1 month at -20 °C.
Associated Products	Anti-hCoV 3504 Anti-hCoV 3505
Buffer	Lyophilized from 50 mM Tris-HCl, pH 7.5, 500 mM NaCl, 0.1 % CHAPS; 6 % sucrose as a stabilizer
Molecular Weight	Predicted molecular weight: 49 kDa
Shelf Life	Unspecified for lyophilized product