

Product Datasheet

Recombinant FABP3, 100 µg 610043

Product Name	Recombinant FABP3, 100 µg
Catalog Number	610043
Description	Recombinant human fatty acid-binding protein 3 (FABP3) with a C-terminal histidine tag.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Storage	Storage at 2-8 °C. After reconstitution stable for 1 month at 2-8°C and 1 month at -20°C.
Associated Products	Anti-h FABP3 2302 Anti-h FABP3 2303 Anti-h FABP3 2304
Buffer	Lyophilized from 50 mM Tris-HCl, pH 7.5; 150 mM NaCl; containing 3 % sucrose, 2 % D-mannitol and 0.01 % Tween 20 as stabilizers.
Molecular Weight	Predicted molecular weight: 16 kDa
Shelf Life	Unspecified for lyophilized product