

Product Datasheet

Human IL-12p70 ELISpot Kit 856.081PC

Product Name	Human IL-12p70 ELISpot Kit
Catalog Number	856.081PC
Description	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
Alternative Names	IL-12
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes natural human IL-12. Specific to p70 heterodimer.
References	Correale J. and Fiol M., Neurology, 2004; 63(12): 2363 - 2370 - Pubmed link Correale, J. et al., Neurology, 2006; 67(4): 652-9. - Pubmed link