

## Product Datasheet

### SAA Antigen, Recombinant 610070

<b>Product Name</b>	SAA Antigen, Recombinant
<b>Catalog Number</b>	610070
<b>Description</b>	Recombinant human serum amyloid A (SAA) protein with an N-terminal histidine tag.
<b>Species/Host</b>	E. coli
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Storage at 2-8 °C. Shelf life after reconstitution 1 month at 2-8 °C and 1 month at -20 °C.
<b>Associated Products</b>	Anti-h SAA clones 2201/RC2201, 2203, 2205, 2208, 2209, 2211, 2212, and 2213
<b>Buffer</b>	Lyophilized from 50 mM Tris-HCl, pH 7.5, 150 mM NaCl, 0.02 % Tween 20; 6 % sucrose as stabilizer
<b>Molecular Weight</b>	Predicted molecular weight: 14 kD
<b>Expression System</b>	E. coli