

Product Datasheet

Recombinant SAA antigen, 100 µg 610070

Product Name	Recombinant SAA antigen, 100 µg
Catalog Number	610070
Description	Recombinant human serum amyloid A (SAA) protein with an N-terminal histidine tag.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Storage	Storage at +2-8 °C. After reconstitution stable for 2 months at +2-8 °C and 2 months at -20 °C.
Associated Products	Anti-h SAA clones 2201/RC2201, 2203, 2205, 2208, 2209, 2211, 2212, and 2213
Buffer	Lyophilized from 50 mM Tris-HCl, pH 7.5, 150 mM NaCl, 0.02 % Tween 20; 6 % sucrose as stabilizer
Molecular Weight	Predicted molecular weight: 14 kD
Shelf Life	24 months
Expression System	E. coli