

## Product Datasheet

### Recombinant PCT antigen, 100 µg 610080

<b>Product Name</b>	Recombinant PCT antigen, 100 µg
<b>Catalog Number</b>	610080
<b>Description</b>	Recombinant human procalcitonin (PCT) protein with a C-terminal histidine tag.
<b>Species/Host</b>	E. coli
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Storage at +2-8 °C. After reconstitution stable for 1 week at +2-8 °C and 10 months at -20 °C.
<b>Associated Products</b>	Anti-h PCT 4003 Anti-h PCT 4004 Anti-h PCT 4005 Anti-h PCT 4006 Anti-h PCT 4008 Anti-h Calcitonin 4051
<b>Buffer</b>	Lyophilized from 50 mM Tris-HCl, pH 7.5; 150 mM NaCl; containing 6 % sucrose as stabilizer
<b>Molecular Weight</b>	Predicted molecular weight: 14 kDa
<b>Shelf Life</b>	Unspecified for lyophilized product
<b>Expression System</b>	E. coli