

Product Datasheet

Human IL-1 beta ELISpot Set 856.101.010P

Product Name	Human IL-1 beta ELISpot Set
Catalog Number	856.101.010P
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 10 non-sterile plates.
Tested Applications	ELISpot
Alternative Names	IL-1b, IL-1beta
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes natural (pro and mature) human IL-1b