

Product Datasheet

Human IL-5 ELISpot Kit 856.111PC

Product Name	Human IL-5 ELISpot Kit
Catalog Number	856.111PC
Description	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-5