

Product Datasheet

NT-proBNP Antigen, Recombinant 610090

Product Name	NT-proBNP Antigen, Recombinant
Catalog Number	610090
Description	Recombinant human N-terminal pro-B-type natriuretic peptide (NT-proBNP) with an N-terminal histidine tag.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Storage	Storage at 2-8°C. Shelf life after reconstitution 1 month at 2-8°C and 1 month at -20 °C
Associated Products	Anti-h NT-proBNP 1306/RC1306, 1307/RC1307, 1308, 1309, 1310, 1311, and 1312
Buffer	Lyophilized from 50 mM Tris-HCl, pH 7.5; 150 mM NaCl; containing 3 % sucrose, 2 % D-mannitol and 0.01 % Tween 20 as stabilizers.
Molecular Weight	Predicted molecular weight: 11 kDa
Expression System	E. coli